

# EASTWIN

www.eastwin.com



*Very handy tool!*

# EQ-6K

## Mini Centrifuge

EQ-6K, Mini centrifuge is mainly used for separating reagent components, collecting tiny droplets to the bottom of tubes and filtering trace reagents. It is an indispensable tool for PCR experiments, suitable for multiple PCR single tubes and 8-strip microtubes. Palm size, easy to move, especially recommended for personal use!

Small and portable

Dual lid-open protection

Quiet and smooth running

Long Life Motor

Rotors Quick-change Without Tools

EASTWIN SCIENTIFIC EQUIPMENTS INC.

# EQ-6K

## Mini Centrifuge

*Very handy tool!*



### Specifications

Rotor	Round Rotor	Tube Strip Rotor
Rotating Speed	6000rpm $\pm$ 5%	
Centrifugal Force	1606 g (0.2ml)	
	1937g (0.5ml)	1990g
	2230 g (1.5ml)	
	2267g (2ml)	
Rotating Direction	Anticlockwise	
Rotor Capacity	8X2ml/1.5ml	4X8-stripX0.2ml
		4X8-stripX0.2ml
		8 single tube
Tube Adapter	8X0.5ml, 8X0.2ml	
Starting Time	2s	
Stopping Time	$\leq$ 2s	
Motor	AC motor, low noise ball bearing	
Power Input	AC 110-250V 50/60Hz	
Dimension L*W*H	150X160X132 mm	
Net Weight	1120g, including round rotor	
Safety Standards	Complying with IEC 61010-2-020	

## EASTWIN

### EASTWIN SCIENTIFIC EQUIPMENTS INC.

Add: F2,Block B,Building 10,#78 Xinglin Street,Suzhou P.R.China

Tel: (+86)0512-85550818 Fax: (+86)0512-89188841

E-maill: info@eastwin.com Web: http://www.eastwin.com



Enterprises Owned by NKY Medical Group (Stock Code 300109)



(+86) 400-086-6869